



# EC DECLARATION OF CONFORMITY

Manufacturer: **WIDEX A/S**  
Nymøllevej 6  
DK-3540 Lyngø  
Denmark

Brand: **WIDEX**  
Product Family: **MAGNIFY**  
Type of Device: **Hearing Aids**  
GMDN Code: **34671 Behind-the-ear air conduction hearing aid**  
**47169 Receiver-in-canal air-conduction hearing aid**  
**34672 In-the-ear air-conduction hearing aid**  
**41209 Canal air-conduction hearing aid**

Product Identification: See next page<s>

We declare under our sole responsibility that above products are in conformity with the following Directives:

---

## **Council Directive 93/42/EEC as amended by Dir. 2007/47/EC (MDD)**

Conformity assessment procedure: Annex II without chapter II.4 of 93/42/EEC

Notified Body: TÜV SÜD Product Service GmbH, Notified Body No.: NB 0123  
Ridlerstr. 65, 80339 München, Germany

Classification of device: Class IIa (according to Annex IX to Council Directive 93/42/EEC)

The products meet the applicable standards and the basic requirements of the Directive 93/42/EEC Annex I.

Applicable standards are listed in the respective technical file.

---

## **Council Directive 2011/65/EU (RoHS) as amended by Dir. 2017/2102/EC (RoHS2)**

Relevant Harmonized Standards: EN 62321

---

## **Council Directive 2014/53/EU (RED)**

Applied standards:  
(with Bluetooth interface) EN 300 328, EN 301 489-1, EN 301 489-17

Applied Standards  
(with Telecoil) EN 300 422-4

Applied Standards  
(with wireless function) EN 300 330, EN 301 489-1, EN 301 489-3

The products are in compliance with the essential requirement of the directive.





# EC DECLARATION OF CONFORMITY

Product Identification	Type of Device
MBB2 (RC, -, M03, M05, M10, M11, M22, M33, M44, DEMO100, DEMO50)	BTE (Behind The Ear) Hearing Aid
MRR2D (RC-mRIC R D, -, M10, M11, M22, M33, M44, DEMO100) MRB0 (RIC 10, -, M03, M05, M10, M11, M22, M33, M44, DEMO100) MRB2D (T-RIC 312 D, -, M05, M10, M11, M22, M33, M44, DEMO100, DEMO50)	RIC (Receiver In the Canal) Hearing Aid
M-XP (T, -, M03, M05, M10, M11, M22, M33, M44) M-IM (L, R, -, M11, M22, M33, M44) M-IP (L, R, -, M11, M22, M33, M44)	ITE (In The Ear) Hearing Aid including faceplate*
M-CIC (L, R, -, M03, M05, M10, M11, M22, M33, M44) M-CIC TR (L, R, -, M03, M05, M10, M11, M22, M33, M44) M-CIC-M (-, M03, M05, M10, M11, M22, M33, M44) M-CIC-M TR (-, M03, M05, M10, M11, M22, M33, M44)	Canal air-conduction hearing aid

\*A Custom made Faceplate is an integral component of the final Product(s) mentioned on this declaration of conformity. Pres-assemblies or configurable variants where CFPS can be made off, are part of the Design Review release process of the Product.